



# **EXHIBIT B**

**This Page Blank (uspto)**

**TeleMed**



# **System Manager's Guide**



**Version 1.5**

**This Page Blank (uspto)**

# System Manager's Guide

**TeleMed**<sup>TM</sup>

**Version 1.5**

**RLIS, Incorporated**

**This Page Blank (uspto)**

Information in this document is subject to change without notice and does not represent a commitment on the part of RLIS, Incorporated.

The software described in this document is furnished under a license agreement. The software may be used or copied only in accordance with the terms of the agreement. Refer to Appendix A, License Agreement, for more details.

Companies, names, and data used in examples herein are fictitious unless otherwise noted.

No part of this document may be reproduced or transmitted in any form or by any means, electronic or mechanical, for any purpose, without the express written permission of RLIS, Incorporated.

© 1998, 1999 RLIS, Incorporated. All rights reserved.

TeleMed, and the TeleMed design (the TeleMed logo) are registered trademarks of RLIS, Incorporated United States of America and other countries.

Microsoft, MS, MS-DOS, Visual FoxPro, FoxPro, FoxPro design (the Microsoft FoxPro logo) are registered trademarks and Windows, and Windows NT are trademarks of Microsoft Corporation in the United States of America and other countries.

Other products mentioned are trademarks or registered trademarks of their respective companies.

This commercial computer software and documentation are being provided to the government with RESTRICTED RIGHTS. Use, duplication, or disclosure by the Government is subject to restrictions as set forth in subparagraph (C)(2) of FAR 52.277-19 or subdivision (c)(1)(ii) of DFAR 52.227-7013. Refer to Appendix A, License Agreement, for more details.

RLIS, Incorporated  
15600 San Pedro, #203  
San Antonio, Texas 78232

Sales and Marketing - 800-496-7547

Document No.: EM-TM001-SM  
Version: 1.5  
Revision: 1.0  
Printed in the United States of America.

*Created by: Daley Professional Group, Inc.*

**Page Intentionally Left Blank**



## Contents

<b>CHAPTER 1: MANAGING THE TELEMED SYSTEM</b>	<b>7</b>
SECURITY	7
Security Overview	7
Telephone Access Security	7
TeleMed System Application Security	7
TeleMed Data Security	8
DATA	8
Storage and Archive Locations	8
TeleMed Data System Features	8
Data Links To Other Hospital Information Systems	8
PRINTERS & PRINTING	9
Printer Configuration	9
Printing with TeleMed	9
MANAGING THE TELEMED SYSTEM	9
Backup	9
System Startup	9
TeleMed Updates/Upgrades	9
STRATEGIES FOR MAINTAINING TELEMED SYSTEM PERFORMANCE	10
Defragmentation:	10
<b>CHAPTER 2: CONFIGURING THE TELEMED SYSTEM</b>	<b>13</b>
TELEMED TASKS AND THE TYPICAL ROLES THAT PERFORM THEM	13
TELEMED ACCOUNT MANAGEMENT	14
Accessing User Rights Screen	14
Creating/Modifying TeleMed User Accounts	15
Deleting TeleMed User Accounts	16
Tips and Hints: TeleMed Account Management	16
MANAGING THE TELEMED PHONE DIRECTORY	17
Creating a Phone Directory Entry	17
Modifying a Phone Directory Entry	18
Deleting a Phone Directory Entry	18
Tips and Hints: Phone Directory Entries	18
MANAGING CUSTOM PRESCRIPTION MENUS	19
Creating Custom Prescription Menus	19
Modifying Custom Prescription Menus	25
Deleting Custom Prescription Menu Entries	27
Tips and Hints: Prescription Menu Entries	28
MANAGING CUSTOM MEDICATIONS MENUS	29
Creating Custom Medication Menus for Administration in the ED	29

<i>Modifying Custom Medication Menus for Administration in the ED</i> .....	35
<i>Deleting Custom Medications Menu Entries</i> .....	37
<i>Tips and Hints: Medications Menu Entries</i> .....	38
MANAGING CUSTOM PREPHRASED TEXT.....	39
<i>Creating Custom Prephrased Text Items</i> .....	39
<i>Modifying Custom Prephrased Text Items</i> .....	40
<i>Deleting Custom Prephrased Text Items</i> .....	41
<i>Tips and Hints: Prephrased Text</i> .....	41
TELEMED DATA MAINTENANCE.....	42
<i>Accessing the Data Maintenance Screen</i> .....	42
<i>Viewing Errors</i> .....	43
<i>Reindex (all) Data</i> .....	43
TESTING TELEMED CHANGES.....	44
<i>Testing TeleMed with a Simulated Patient</i> .....	44
<b>CHAPTER 3: TELEMED SUPPORT SERVICES</b> .....	<b>45</b>
TELEMED™ SERVICE.....	45
CONTACTING TELEMED HELP DESK.....	46
<b>APPENDIXES</b> .....	<b>47</b>
APPENDIX A: TELEMED LICENSE AGREEMENT.....	47
<b>INDEX</b> .....	<b>51</b>

## Chapter 1: Managing the TeleMed System

### Security

#### Security Overview

In order to maintain TeleMed system's security protocol, prudence dictates that we not publicize how security is maintained. Suffice it to say, your TeleMed system has utilized a combination of security measures designed to meet your individual operation's security needs. These measures include some or all of the following:

- Modem security features of NT Server (callback, password, encryption or silent modems)
- Commercial firewall products
- Password schemes (multiple very long passwords, frequently automatically changing passwords, one-time-use passwords, etc.)
- NT's system level security (up to the C2 level of security required by many government agencies)

Your TeleMed installer coordinated with your hospital's MIS representative on actual security configuration for your TeleMed system. Upon completion of the TeleMed installation, your operation was provided a TeleMed Installation Specification that describes how system security was implemented.

#### Telephone Access Security

The TeleMed system can accept transcription file returns and maintenance access by phone. The only function allowed by TeleMed during transcription access is the deposit (upload) of transcription files.

Maintenance/support dial-in access is controlled by a combination of operating system security features and firewall products which prevent unauthorized access by using sophisticated techniques that include some or all of industry standard security schemes previously listed.

#### TeleMed System Application Security

Access to the TeleMed system application is restricted to individuals whose names have been entered, in advance, to the list of eligible active users. Each of these users must also enter a private password, which grants access to components of the software appropriate for that user. (For example, only the physician may edit his or her own portion of the record. Nursing notes are similarly proprietary.) After a certain period of inactivity, terminals revert to a "locked" status which requires re-entry of passwords, etc. in order to function.

## **TeleMed Data Security**

[TBS]

## **Data**

### **Storage and Archive Locations**

Upon completion of the TeleMed installation, your operation was provided a TeleMed Installation Specification that describes where your TeleMed system's data is stored and archived (backed-up).

### **TeleMed Data System Features**

- The TeleMed system is built on an X-Base data structure which allows other "off the shelf" database products to query and create other "operation specific" reports.
- The TeleMed system operates on a Windows NT operating system (version 4.0 or higher) thus providing an industry standard operating system and network platform.
- The TeleMed system is Network Operating System (NOS) independent allowing flexibility for future information system iterations.
- The TeleMed system supports Structured Query Language (SQL) based query of patient medical data supporting industry standard data exchange interface.
- The TeleMed system utilizes "Rushmore" search optimization technology to meet system response demands of a Level I Trauma Center.

### **Data Links To Other Hospital Information Systems**

Your TeleMed installer coordinated with your hospital's MIS representative on actual security configuration for your TeleMed system. Upon completion of the TeleMed installation, your operation was provided a TeleMed Installation Specification that describes how system security was implemented.

## **Printers & Printing**

### **Printer Configuration**

TeleMed is a Windows (NT 4.0 or 95) application and conforms to Windows NT 4.0 client-server printer sharing features. Your TeleMed system printer(s) can be configured to operate as directly connected network interface print devices or via a print server.

Windows NT 4.0 Print Manager takes care of all printer maintenance (at an operating system level).

Your TeleMed installer tailored your printer configuration to meet your ED and system management needs. Upon completion of the TeleMed installation, your operation was provided a *TeleMed Installation Specification* that describes how printer management was implemented.

### **Printing with TeleMed**

TeleMed provides its own printer interface. This interface allows for a variety of record and report printing. Refer to your *TeleMed User's Guide* index (*Printing*) for instructions on how to access these print functions.

## **Managing The TeleMed System**

Management demands of TeleMed Servers differs for each installation. This section provides a general description of common system management activities.

### **Backup**

[TBS]

### **System Startup**

Typically, TeleMed is configured to startup automatically upon NT Workstation and Server restart. Some personalized workstations may not have TeleMed configured to auto start.

### **TeleMed Updates/Upgrades**

RLIS distributes quarterly updates of TeleMed to their Warranty and Support customers. These distributions include

- **Maintenance Releases** are minor software enhancements.
- **Quarterly Updates** are provided to implement changes in prescription, standard forms or medical information. These updates will be distributed in the following months: January, April, July, October. Quarterly Updates must be installed within 45 days of distribution.

**WARNING:** Failure to implement Quarterly Updates within the 45 day time period will result in automatic shutdown of the prescription and medications components of the TeleMed system.

**NOTE:** In the interest of practicing good medicine, these quarterly updates address: pharmaceutical database updates, 3<sup>rd</sup> party licensing terms, medical knowledge base updates and other enhancements as required for prudent medical practice.

- **Periodic Upgrades** are distinctly new versions of TeleMed and are included in all TeleMed Service options (Platinum, Gold and Silver). These upgrades update all the TeleMed licensed products currently installed on your TeleMed system.
- **Functional Upgrades** are upgrades for separate optional TeleMed functions. These upgrades are sold as separately licensed TeleMed modules. Once purchased and installed, these functional upgrades are updated through periodic upgrades.

Each quarterly update will contain installation instructions.

Your TeleMed installer tailored your TeleMed system to meet your ED and system management needs. Upon completion of the TeleMed installation, your operation was provided a TeleMed Installation Specification that describes a recommended TeleMed/NT system management schedule.

## ***Strategies For Maintaining TeleMed System Performance***

All TeleMed systems are installed with a suite of Windows NT specific system/disk utilities geared toward optimizing operating system performance and reliability. This section is not intended to be a "how to" on NT performance optimizations. We do wish to emphasize the importance of the background defragmentation utility deployed with TeleMed system disk drives and files. The following section (Defragmentation) contains excerpts from a white paper<sup>1</sup> that describes the effects of defragmentation on NT system performance.

### **Defragmentation:**

**Background:** The Windows NT File System (NTFS) was created specifically for use with the Windows NT Operating System. NTFS is a recoverable, secure, reliable and efficient file system which supports client-server systems.

NTFS keeps track of the contents of a volume in a file called the Master File Table (MFT). The MFT consists of an array of file records, with each file usually represented in the MFT by one file record. If, however, a file has a large number of file attributes, such as security information, or becomes highly fragmented, more than one file record may be needed. As more files and file records are added to a volume, the MFT may need to expand to accommodate more file records. The ability of the MFT to expand is part of the way NTFS supports large disks with large numbers of files. As the MFT expands, it may also fragment across the disk. It should be noted that once the MFT expands it will not shrink, but as files are deleted from the volume, unused file records become available for use by NTFS as new files are subsequently added.

---

<sup>1</sup> authored by Executive Software (makers of Diskeeper). The complete report can be found on Executive Software's web site ([www.execsoft.com](http://www.execsoft.com)).

The tests clearly show that performance of NTFS suffers from fragmentation of files and free space. The degree of performance degradation is substantial when the volume is mildly fragmented and worsens steadily as fragmentation increases. As can be seen in all the tests, defragmentation of files and free space significantly restores system performance. In fact, defragmenting the volume almost completely restores disk performance to pre-fragmentation levels.

Comparison of the defragmented volume to the control volume does show, however, that even completely defragmenting files cannot fully restore the performance lost because of fragmentation. This is because of the fragmentation of the Master File Table. The amount of time to read and write files to the volume did not return to pre-fragmentation levels, even when transfer rates did. This shows that a fragmented Master File Table requires additional reads from the volume and shows that keeping a disk volume defragmented from the time it is initialized is very important. Failure to keep an NTFS disk volume defragmented will result in progressively slower file operations. If allowed to continue so long that the Master File Table itself becomes fragmented, this can result in a volume that cannot be restored to full performance even with complete defragmentation.

Testing indicates that fragmentation is a major factor in the performance of a Windows NT system. Defragmentation of an NTFS volume, especially one that has been in use for some length of time, can be expected to substantially improve the overall responsiveness of the system.

Your TeleMed installer tailored your Windows NT system/disk utilities to meet your ED and system management needs. Upon completion of the TeleMed installation, your operation was provided a TeleMed Installation Specification that describes how printer management was implemented.





## Chapter 2: Configuring the TeleMed System

### ***TeleMed Tasks And The Typical Roles That Perform Them***

The following chart lists some of the typical roles and associated tasks they can perform with TeleMed.

	Doctor	Nurse	Tech	Clerk	Comment:
Enter New Patient	•	•			
View Clinical Summaries	•	•			
Edit Doctor Drug Preferences	•				
Create Prescriptions	•				
Order Medications	•	•			
Order Therapeutics	•	•			
Generate Reports	•	•		•	
Quit Program	*				Optional
Edit Users	*				Optional
System Manager	*	*			Optional
Delete Insurance Plans	*				Optional
Edit Prephrased Text	*				Optional
Sign Physician Med. Record	*				Optional
Sign Nurse Medical Record		*			Optional


**Table 1: Typical TeleMed Tasks & Roles Chart**

**NOTE:** The tasks that an individual is authorized to perform in TeleMed are driven by your medical operation's policies and procedures. This chart is intended to familiarize those new to TeleMed and is not intending to influence your medical operation's policies or procedures.

## TeleMed Account Management

### Accessing User Rights Screen

From the "Active Patient List" screen



1. Click the Utilities button.
2. Click on the Edit Users button.
3. Enter your password in the password prompt window
4. Click on the Continue button.

**User Rights**

Full Name (First, MI, Last) (30 Characters Max)

Short Name For Button Selection (20 Characters Max)

Doctor State Number

Doctor Federal Number

Doctor Telephone Number

Password (8 Characters Max)

PIN (8 Characters Max)

☒ Enter New Patient  
☒ View Clinical Summaries  
☒ Edit Doctor Drug Prefs  
☒ Create Prescriptions  
☒ Order Medications  
☒ Order Therapeutics  
☒ Generate Reports  
☒ Quit Program  
☒ Edit Users  
☒ System Manager

☒ Delete Insurance Plans  
☒ Edit Prephased Test  
☐ Sign Physician Med Record  
☐ Sign Nurse Medical Record

Previous  
 Next  
 Go To  
 Change  
 Delete  
 Exit

To Select A Name or Position, Double Click on Selection In Button Position Column  
 or  
 Click on Selection Then Click "Go To"

Button Position	Name or Position
	Olivares, UC
	LaPoint, CNA II
	Turner, CNA II
	Campos, CNA II
	Esparza, CNA
	Perez, CNA II
	Registration 1
	Daley, RN
	Garza, RN
	Scott, RN
	13
	14
	Hawkins, RN
	Moreno, LVN
	Rodriguez, LVN
	J. Duque, RN
	Smith, RN
	Ross
	Townsend, RN
	Bradshaw, RN
	Wilkins, RN
	Flores, RN
	Phillips, RN
	26
	Orgill, RN
	Cynthia, RN
	Sanders, RN
	Cervantes, RN
	31
	Youngblood, RN
	Carroll, RN

**Figure 1: User Rights**

(Image: user-rights-er0157.bmp)

### Accessing User Rights Screen (continued)

**Note:**

- On the right is a list of all current users.
- The **Full Name** is the name displayed in most places and printed.
- The **Short Name** is the name displayed on the user list and the log on buttons.
- The **Password** is limited to eight characters and is not case sensitive.
- The boxes in the lower left are user rights.
  - \* Selecting **System Manager** automatically checks several boxes (as shown in Figure 1: User Rights and Figure 2: Default System Manager Rights).
- After additions or changes are made, the Change button must be selected to store the changes.

### Creating/Modifying TeleMed User Accounts


#### *From the "User Rights screen"*



1. Select the account name from the user list at the right of screen
2. Click the Go To button. (or Double-click in account name)
3. Click on check boxes and modify form fields as appropriate.
4. After entries are made, click the Change button must be selected to store the changes.
5. Click the Exit button to return to the Utilities screen
6. Click on the Return to Program button to log back into TeleMed.

## Deleting TeleMed User Accounts

*From the "User Rights screen"*



1. Select the account name from the user list at the right of screen. You can use the Previous / Next buttons or click directly on the account name.
2. Click the Go To button. (or Double-click in account name)
3. Click on Delete button. This will delete (blank out) all fields and check boxes except the **Short Name**. The **Short Name** entry field will revert to a default entry number with no rights. Modify this null account as appropriate
4. Click the Exit button to return to the Utilities screen
5. Click on the Return to Program button to log back into TeleMed.

## Tips and Hints: TeleMed Account Management

- System Manager privileges are assigned to a user much like one is granted privileges to "Generate Reports" (refer to Table 1: Typical TeleMed Tasks & Roles Chart).
- When the "System Manager" rights box is selected several boxes are automatically checked (selected):

- |  |  |
|--|--|
| <input checked="" type="checkbox"/> Enter New Patient            | <input checked="" type="checkbox"/> Delete Insurance Plans |
| <input checked="" type="checkbox"/> View Clinical Summaries      | <input checked="" type="checkbox"/> Edit Prephrased Text   |
| <input checked="" type="checkbox"/> Edit Doctor Drug Preferences | <input type="checkbox"/> Sign Physician Med. Record        |
| <input checked="" type="checkbox"/> Create Prescriptions         | <input type="checkbox"/> Sign Nurse Medical Record         |
| <input checked="" type="checkbox"/> Order Medications            |  |
| <input checked="" type="checkbox"/> Order Therapeutics           |  |
| <input checked="" type="checkbox"/> Generate Reports             |  |
| <input checked="" type="checkbox"/> Quit Program                 |  |
| <input checked="" type="checkbox"/> Edit Users                   |  |
| <input checked="" type="checkbox"/> System Manager               |  |

**Figure 2: Default System Manager Rights**

## Managing the TeleMed Phone Directory

TeleMed's "Phone Directory" function allow you to quickly enter, modify and retrieve any referral, consultant, service or other medical organization from a common interface.

### Creating a Phone Directory Entry

*From the "Active Patient List" or "Medical Information" screen*

**Phone Directory**

**New - Add**

1. Click the Phone Directory button at the right of screen.
2. Click on the New - Add button.

**Personnel Directory** Active User: Marcus Welby

Name: Last  First  MI  Generation

Firm

Address

City  State  Zip

Office Phone  Service Phone

Alt. Phone  Home Phone

Beeper  Car Phone

Fax

☒ Digital Pager ☐ Voice Pager

Medical Certification  ☐ Notify Before Treatment Notify From  To

Medical Specialty

Managed Care Systems

Referral Preferences

Phone	Doctor	Specialty
(210) 377-3366	Dougie Houser	Pediatrics
(210) 377-3366	Bebe Newbirth	OB-GYN
(210) 377-3366	Bones T McCoy	Orthopedics
(210) 377-3366	Ben Casey	General Surgery
(210) 377-3366	Arthur C Hart	Cardiology

**On Call**

Specialty

Start Date

Start Time

End Date

End Time

Time must be 24-hour format

Clear

Save New Entry

Save Changes

Delete

List Directory

Clear User

Exit

**Figure 3: Personnel Directory**

(Image: personnel-directory-cr0132.bmp)

**Save New Entry**


**Save Changes**

or

3. Click on check boxes and complete form fields as appropriate.
4. Click on the Save New Entry (end then click on the Exit button)  
or just click on the Save Changes button at the right of screen.

### **Modifying a Phone Directory Entry**


*From the "Active Patient List" or "Medical Information" screen*



1. Click the Phone Directory button at the right of screen.
2. Click on the Show Detail button. This will open the Personnel Directory screen (refer to Figure 3: Personnel Directory, page 17).
3. Click on check boxes and modify form fields as appropriate.
4. Click on the Save Changes button at the right of screen.

### **Deleting a Phone Directory Entry**

*From the "Active Patient List" or "Medical Information" screen*



1. Click the Phone Directory button at the right of screen.
2. Click on the Show Detail button. This will open the Personnel Directory screen (refer to Figure 3: Personnel Directory, page 17).
3. Click on the Delete button at the right of screen. (the Delete Entry warning window will pop up)
4. Click on the Delete button in the Delete Entry warning window.

### **Tips and Hints: Phone Directory Entries**

- If you ever need to enter two separate directory entries for one person (or organization), use the Save New Entry button. This will allow you to maintain reuse any common information without retyping.

## Managing Custom Prescription Menus

### Creating Custom Prescription Menus

*From the "Active Patient List" screen*



1. Click the Utilities button.




2. Click on the Edit Doctor Drug Prescription Pref. List button.



3. Enter your password in the password prompt window

4. Click on the Continue button.



**Figure 4: Custom Drugs (prescription menu)**

(Image: custom-drugs\_prescription-cr0134.bmp)



5. Click the Add button.

**Creating Custom Prescription Menus (continued)****Option 1: By Drug Name**

6. Enter the first two (or more) letters of the drug name.



7. Click on the Find button.

**Drug Alphabetical**

Active User Marcus G. Welby.

Patient Name	Adm #	Pat #	SS #	Bed
<div><div><div>Zitamin (Prenatal Multivitamins &amp; Minerals w/ Iron &amp; FA Tab 1 MG)</div><div>Zithromax 250 MG (Azithromycin Cap 250 MG)</div></div><div></div></div> <div><div>Enter First 2 or more Characters of Drug Name</div><div>ZIT</div><div>Find</div><div>Classification Drug List</div><div>2 2</div><div>Clear User</div><div>Exit</div></div>				

**Figure 5: Drug Alphabetical**

(Image: drug-alphabetical\_zit-er0139.bmp)

8. Select (*Double click*) the drug.



**Creating Custom Prescription Menus (continued)****Option 1: By Drug Name (continued)**

No.	Type	Route	Frequency	Duration	Name Displayed on Your Custom List
1/2	tab	po	1 time only	1 day	Zithromax 250 MG
3/4	cap	pr	q day	2 days	Drug Name
1	inck	os	bid	3 days	Zithromax 250 MG
2	lozenge	od	tid	4 days	Generic
3	anpuke	ou	qid	5 days	Azithromycin Cap 250 MG
4	packet	as	5 times daily	7 days	Drug Group
5	sup	ad	q 3'	10 days	Drug Subgroup
6	preco	au	q 3-4'	12 days	Drug Class
7	implant	inhalations	q 3-6'	14 days	
8	patch	intranasal	q 3-12'	21 days	
9	bar	apply to affec. area	q 4'	28 days	
10	bottle	topically	q 4-6'	30 days	
11	qtts	sublingual	q 4-8'		Quantity 10
+ 1/2	tsp	vaginal	q 5'		Duration 5 Days Remaining 123
lbs	as ingestion	q 6'			SIG (Use Option Return to manually break lines)
cc	transdermal	q 6-8'			2 caps today, then 1 cap daily on empty stomach
mcg	IM	q 6-12'			
mg	IV	q 8'			
ga	subcutaneous	q 8-12'			
ml		q 12'			
mlu					
mu					
iu					
units					

☐ + 1/2

☐ As Needed

To add 1/2 to the No. (i.e. 3-1/2) click the box on "+ 1/2" then click the appropriate No. button.

SIG may be built using buttons or typed.

The same drug can be repeated with different specifications. Use a different name for your custom list.

Save Change  
Save New  
Clear User  
Exit

**Figure 6: Zithromax 250 Drug Screen**

(Image: zithromax\_250-er0149.bmp)

9. Edit the display name, drug name, generic, quantity, duration and SIG fields as required.

Save New

10. Click the Save New. (returns to custom drug list)

Exit

11. Click the Exit button.

## Creating Custom Prescription Menus (continued)

### Option 2: By Drug Classification

Classification Drug List

6. Click on the Classification Drug List button.

**Drug Groups**

Analgesic & Anesthetics	Biologicals	Endocrine & Metabolic Drugs	Hematological Agents	Misc. Products
Analgesics Narc.	Antisera	Androgen-Anabolic	Anticoagulants	Antidotes
Analgesics NonNarc.	Toxoids	Antidiabetic	Hematopoietic	Antiseptic & Disinfect
Antirheumatic	Vaccines	Contraceptives	Hemostatics	Chemicals
General Anesthetics	Misc. Biologicals	Corticosteroids	Misc. Hematological	Diagnostic Products
Gout		Estrogens	Neuromuscular	Medical Devices
Local Anesthetics	Cardiovascular Agents	Oxytocics	Anticonvulsant	Pharm. Adjuvants
Migraine Products	Antianginal Agents	Progestins	Antimaysthenic	Unclassified
Antinfective Agents	Antiarrhythmic	Thyroid	Antiparkinsonian	
Anoebicides	Antihyperlipidemic	Misc. Endocrine	Neuromusc. Blockers	
Aminoglycosides	Antihypertensive		Skel. Musc. Relaxants	
Anthelmintic	Beta Blockers	GI Agents		Nutritional Products
Antifungals	Calcium Blockers	Antacids		Dietary Products
Antimalarial	Cardiotonics	Antidiarrheals	Respiratory Agents	Mineral & Electrolyte
Antitubercular	Diuretics	Antiemetics	Antisthmatics	Multivitamins
Antiviral	Pressors	Digestive Aids	Antihistamines	Nutrients
Cephalosporins	Misc. Cardiovascular	Laxatives	Decongestants	Vitamins
Fluoroquinolones		Ulcer Drugs	Cough/Cold	
Macrolide Antibiotics	Central Nervous Sys.	Misc. GI	Misc. Respiratory	
Penicillins	Antianxiety Agents	GU Products		
Sulfonamides	Antidepressants	Urinary Antinfectives	Topical Products	
Tetracyclines	Antipsychotics	Urinary Antispasmodics	Dermatological	
Misc. Antinfectives	Hypnotics	Vaginal Products	Mouth & Throat	
Antineoplastic Agents	Stimulants	Misc. GU Products	Ophthalmic	
Antineoplastics	Misc. Psychother.		Otic	
				Clear User
				Exit

Figure 7: Drug Class

(Image: drug-groups-cr0141.bmp)

7. Click on the appropriate drug group button.  
(ex. Macrolide Antibiotics)

**Creating Custom Prescription Menus (continued)****Option 2: By Drug Classification (continued)**

**Drug Picker** Active User: Thomas Daley

Patient Name:  Adm #  Pat #  SS #  Bed

Clindamycin HCl 150 MG (Clindamycin HCl Cap 150 MG)  
 Cleocin 300 MG (Clindamycin HCl Cap 300 MG)  
 Clindamycin HCl 300 MG (Clindamycin HCl Cap 300 MG)  
 Clindamycin HCl (Clindamycin HCl Powder)  
 Cleocin 75 M G/5ML (Clindamycin Palmitate For Susp 75 MG/5ML)  
 Clindamycin Palmitate 75 MG/5ML (Clindamycin Palmitate For Susp 75 M)  
 Cleocin 150 MG/ML (Clindamycin Phosphate Inj 150 MG/ML)  
 Clindamycin Phosphate 150 MG/ML (Clindamycin Phosphate Inj 150 MG/  
 Clindamycin Phosphate (Clindamycin Phosphate Powder)  
 Cleocin (300 ) (Clindamycin Phosphate in D5W IV Soln 300 MG/50ML)  
 Clindamycin Phosphate in D5W 300 MG/50ML (Clindamycin Phosphate in  
 Cleocin (600 ) (Clindamycin Phosphate in D5W IV Soln 600 MG/50ML)  
 Clindamycin Phosphate in D5W 600 MG/50ML (Clindamycin Phosphate in  
 Cleocin (900 ) (Clindamycin Phosphate in D5W IV Soln 900 MG/50ML)  
 Clindamycin Phosphate in D5W 900 MG/50ML (Clindamycin Phosphate in  
 Azithromycin 250 MG (Azithromycin Cap 250 MG)  
**Zithromax 250 MG (Azithromycin Cap 250 MG)**  
 Biaxin 250 M G (Clarithromycin Tab 250 MG)  
 Clarithromycin 250 MG (Clarithromycin Tab 250 MG)  
 Biaxin 500 M G (Clarithromycin Tab 500 MG)  
 Clarithromycin 500 MG (Clarithromycin Tab 500 MG)  
 Biaxin 125 M G/5ML (Clarithromycin For Susp 125 MG/5ML)  
 Clarithromycin 125 MG/5ML (Clarithromycin For Susp 125 MG/5ML)  
 Biaxin 250 M G/5ML (Clarithromycin For Susp 250 MG/5ML)  
 Clarithromycin 250 MG/5ML (Clarithromycin For Susp 250 MG/5ML)

**Drug Group**  
 Antimicrobial Agents

**Drug Subgroup**  
 Macrolide Antibiotics

**Drug Class**  
 03 109

**Figure 8: Drug Picker**

(Image: drug-picker\_zit-er0143.bmp)

8. Select (*Double click*) the drug.

**Creating Custom Prescription Menus (continued)****Option 2: By Drug Classification (continued)**

No.	Type	Route	Frequency	Duration	Name Displayed on Your Custom List
<input checked="" type="checkbox"/> 1/4	tab	po	1 time only	1 day	Zithromax 250 MG
<input type="checkbox"/> 1/2	cap	pr	q day	2 days	Drug Name
<input type="checkbox"/> 3/4	inch	as	bid	3 days	Zithromax 250 MG
<input type="checkbox"/> 1	lozenge	od	tid	4 days	Generic
<input type="checkbox"/> 2	ampule	ou	qid	5 days	Azithromycin Cap 250 MG
<input type="checkbox"/> 3	packet	as	5 times daily	7 days	Drug Group
<input type="checkbox"/> 4	sup	ad	q 3*	10 days	Drug Subgroup
<input type="checkbox"/> 5	piece	ar	q 3-4*	12 days	Drug Class
<input type="checkbox"/> 6	inhalant	inhalations	q 3-6*	14 days	Quantity <input type="text" value="10"/>
<input type="checkbox"/> 7	patch	intranasal	q 3-12*	21 days	Duration <input type="text" value="5"/> Days Remaining <b>123</b>
<input type="checkbox"/> 8	bar	apply to affec. area	q 4*	28 days	SIG (Use Option Return to manually break lines)
<input type="checkbox"/> 9	bottle	topically	q 4-6*	30 days	2 caps today, then 1 cap daily on empty stomach
<input type="checkbox"/> 10	qts	sublingual	q 4-8*		
<input type="checkbox"/> + 1/2	lsp	vaginal	q 5*		
	tbz	as irrigation	q 6*		
	cc	transdermal	q 6-8*		
	mcg	IM	q 6-12*		
	mg	IV	q 8*		
	gm	subcutaneous	q 8-12*		
	ml		q 12*		
	mcu				
	mu				
	ru				
	units				

☐ As Needed

To add 1/2 to the No. (i.e. 3-1/2) click the box on "+ 1/2" then click the appropriate No. button.

SIG may be built using buttons or typed.

The same drug can be repeated with different specifications N use a different name for your custom list.

Save Change  
Save New  
Clear User  
Exit

**Figure 9: Zithromax 250 Drug Screen**

(Image: zithromax\_250-er0149.bmp)

9. Edit the display name, drug name, generic, quantity, duration and SIG fields as required.

**Save New**

10. Click the Save New. (returns to custom drug list)

**Exit**

11. Click the Exit button.

## Modifying Custom Prescription Menus

From the "Active Patient List" screen



1. Click the Utilities button.



2. Click on the Edit Doctor Drug Prescription Pref. List button.



3. Enter your password in the password prompt window  
4. Click on the Continue button.




Figure 10: Custom Drugs (prescription menu)

(Image: custom-drugs\_prescription-er0134.bmp)

5. Select (*click on*) the drug to be modified from the drug list at left of screen.



6. Click on the Change button.

### Modifying Custom Prescription Menus (continued)

No.	Type	Route	Frequency	Duration	Name Displayed on Your Custom List
<input checked="" type="checkbox"/> 1/4	tab	po	1 time only	1 day	Zithromax 250 MG
<input type="checkbox"/> 1/2	cap	pr	q day	2 days	Drug Name
<input type="checkbox"/> 3/4	inch	os	bid	3 days	Zithromax 250 MG
<input type="checkbox"/> 1	lozenge	od	tid	4 days	Generic
<input type="checkbox"/> 2	ampule	ou	qid	5 days	Azithromycin Cap 250 MG
<input type="checkbox"/> 3	packet	as	5 times daily	7 days	Drug Group
<input type="checkbox"/> 4	sup	ad	q 3*	10 days	Drug Subgroup
<input type="checkbox"/> 5	piece	au	q 3-4*	12 days	Drug Class
<input type="checkbox"/> 6	implant	inhalations	q 3-6*	14 days	Quantity <input type="text" value="10"/>
<input type="checkbox"/> 7	patch	intranasal	q 3-12*	21 days	Duration <input type="text" value="5"/> Days Remaining <b>123</b>
<input type="checkbox"/> 8	bar	apply to affec. area	q 4*	28 days	SIG (Use Option Return to manually break lines)
<input type="checkbox"/> 9	bottle	topically	q 4-6*	30 days	2 caps today, then 1 cap daily on empty stomach
<input type="checkbox"/> 10	qtts	sublingual	q 4-8*		
<input type="checkbox"/> + 1/2	tsp	vaginal	q 5*		
	lbs	as ingestion	q 6*		
	cc	transdermal	q 6-8*		
	mcg	IM	q 6-12*		
	mg	IV	q 8*		
	qn	subcutaneous	q 8-12*		
	ml		q 12*		
	miu			<input type="checkbox"/> As Needed	
	mu				
	iu				
	units				

To add 1/2 to the No. (i.e. 3-1/2) click the box on "+ 1/2" then click the appropriate No. button.

SIG may be built using buttons or typed.

The same drug can be repeated with different specifications. Use a different name for your custom list.

**Figure 11: Zithromax 250 Drug Screen**

(Image: zithromax\_250-er0149.bmp)


- Edit the display name, drug name, generic, quantity, duration and SIG fields as required.
- Click the Save Change button. (returns to custom drug list)
- Click the Exit button.

## Deleting Custom Prescription Menu Entries

*From the "Active Patient List" screen*



1. Click the Utilities button.




2. Click on the Edit Doctor Drug Prescription Pref. List button.



3. Enter your password in the password prompt window

4. Click on the Continue button.



**Figure 12: Custom Drugs (prescription menu)**


(Image: custom-drugs\_prescription-er0134.bmp)

5. Select (*click on*) the drug to be deleted from the drug list at left of screen.

6. Click on the Delete button.



7. Click on the Exit button.



### **Tips and Hints: Prescription Menu Entries**

- A different version of a drug can also be entered and Save New will save the additional version without changing the previous version.
- Note that 165 drugs can be listed in the custom drug list. The Remaining Available Positions notes how many positions remain.



## Managing Custom Medications Menus

### Creating Custom Medication Menus for Administration in the ED

*From the "Active Patient List" screen*

Utilities

1. Click the Utilities button.

Edit Doctor Drug Order Pref. List

2. Click on the Edit Doctor Drug Order Pref. List button.

Continue

3. Enter your password in the password prompt window

4. Click on the Continue button.

**Custom Drugs** Active User: Marcus G. Welby, MD

Alupent Inhaler, 2 puffs  
Biaxin 500 MG  
Cefclor 250 MG  
Cefclor 250 MG/5ML  
DT Adult  
DT Pediatric dose, 6.6-5 LFU  
Demerol 50 MG IM  
Doxycycline 100 MG  
Duricef 500 MG  
Flagyl 2 gm for trich  
Naprosyn 500 MG  
Nubain 10 MG IM  
Orudis 75 MG  
Penicillin G Pot Inj 1 Mu/D5W  
Phenergan 12.5 MG IV  
Phenergan 25 MG IM  
Phenergan Supp 25 MG  
Proventil Inhaler, 2 puffs  
Proventil by nebulizer  
Relafen 500 MG  
Rocephin 1 gm IM  
Rocephin 1 gm IV  
Rocephin 250 mg IM  
Td, adult  
Tylenol #3, 1 tab po  
Tylenol #3, 2 tabs po

Add  
Delete  
Change

164 Remaining Available Positions

Clear User  
Exit

**Figure 13: Custom Drugs (drug order menu)**

(Image: custom-drugs\_order-er0134.bmp)


Add

5. Click the Add button.

## Creating Custom Medication Menus for Administration in the ED (continued)

### Option 1: By Drug Name

6. Enter the first two (or more) letters of the drug name.



7. Click on the Find button.

### Drug Alphabetical

Active User Marcus G. Welby.

Patient Name

Adm #

Pat #

SS #

Bed

A screenshot of a software window titled "Drug Alphabetical". The window is divided into two main sections. The left section is a large list box containing two items: "Zitatin [Prenatal Multivitamins & Minerals w/ Iron & FA Tab 1 MG]" and "Zithromax 250 MG (Azithromycin Cap 250 MG)". The right section contains a search interface. At the top, it says "Enter First 2 or more Characters of Drug Name" above a text input field containing "ZIT". Below the input field is a "Find" button. Underneath the "Find" button is a "Classification Drug List" button. At the bottom right of the window are two buttons: "Clear User" and "Exit".

**Figure 14: Drug Alphabetical**

(Image: drug-alphabetical\_zit-er0139.bmp)

8. Select (*Double click*) the drug.

## Creating Custom Medication Menus for Administration in the ED (continued)

### Option 1: By Drug Name (continued)

No.	Type	Route	Frequency	Duration	Name Displayed on Your Custom List
1/4	tab	po	1 time only	1 day	Zithromax 500 mg load dose
1/2	cap	pr	q day	2 days	Drug Name
3/4	mch	os	bid	3 days	Zithromax 250 MG
1	lozenge	od	tid	4 days	Generic
2	ampule	ou	qid	5 days	Azithromycin Cap 250 MG
3	packet	as	5 times daily	7 days	Drug Group
4	sip	ad	q 3"	10 days	Drug Subgroup
5	piece	au	q 3-4"	12 days	Drug Class
6	explant	inhalations	q 3-6"	14 days	
7	patch	intranasal	q 3-12"	21 days	
8	bar	apply to affect. area	q 4"	28 days	
9	bottle	topically	q 4-6"	30 days	
10	qtts	sublingual	q 4-8"		Quantity <input type="text" value="10"/>
<input type="checkbox"/> + 1/2	tsp	vaginal	q 5"		Duration <input type="text" value="5"/> Days Remaining <input type="text" value="165"/>
	tbs	as irrigation	q 6"		SIG (Use Option Return to manually break lines)
	cc	transdermal	q 6-8"		<input type="text" value="2 cap po now"/>
	mcc	IM	q 6-12"		
	mg	IV	q 8"		
	gm	subcutaneous	q 8-12"		
	ml		q 12"		
	mu				
	au				
	unit				

☐ As Needed

To add 1/2 to the No. (i.e. 3-1/2) click the box on "+ 1/2" then click the appropriate No. button.

SIG may be built using buttons or typed.

The same drug can be repeated with different specifications. Use a different name for your custom list.

Save Change  
Save New  
Clear User  
Exit

**Figure 15: Zithromax 500 mg Drug Screen**

(Image: zithromax\_500-er0149.bmp)

9. Edit the display name, drug name, generic, quantity, duration and SIG fields as required.

Save New

10. Click the Save New. (returns to custom drug list)

Exit

11. Click the Exit button.

## Creating Custom Medication Menus for Administration in the ED (continued)

### Option 2: By Drug Classification

**Classification Drug List**

6. Click on the Classification Drug List button.

**Drug Groups**

Analgesic & Anesthetics	Biologicals	Endocrine & Metabolic Drugs	Hematological Agents	Misc. Products
Analgesics-Narc.	Antisera	Androgen Anabolic	Anticoagulants	Antidotes
Analgesics-NonNarc.	Toxoids	Antidiabetic	Hematopoietic	Antiseptic & Disinfect
Antirheumatic	Vaccines	Contraceptives	Hemostatics	Chemicals
General Anesthetics	Misc. Biologicals	Corticosteroids	Misc. Hematological	Diagnostic Products
Gout		Estrogens	Neuromuscular	Medical Devices
Local Anesthetics	Cardiovascular Agents	Oxytocics	Anticonvulsant	Pharm. Adjuvants
Migraine Products	Antianginal Agents	Progestins	Antiparkinsonian	Unclassified
Antianfective Agents	Antiarrhythmic	Thyroid	Neuromusc. Blockers	
Amoebicides	Antihyperlipidemic	Misc. Endocrine	Skel. Musc. Relaxants	
Aminoglycosides	Antihypertensive	GI Agents		Nutritional Products
Anthelmintic	Beta Blockers	Antacids		Dietary Products
Antifungals	Calcium Blockers	Antidiarrheals	Respiratory Agents	Mineral & Electrolyte
Antimalarial	Cardiotonics	Antiemetics	Antisrhetics	Multivitamins
Antitubercular	Diuretics	Digestive Aids	Antihistamines	Nutrients
Antiviral	Pressors	Laxatives	Decongestants	Vitamins
Cephalosporins	Misc. Cardiovascular	Ulcer Drugs	Cough/Cold	
Fluoroquinolones		Misc. GI	Misc. Respiratory	
Macrolide Antibiotics	Central Nervous Sys.			
Penicillins	Antianxiety Agents	GU Products	Topical Products	
Sulfonamides	Antidepressants	Urinary Antinfectives	Dermatological	
Tetracyclines	Antipsychotics	Urinary Antispasmodics	Mouth & Throat	
Misc. Antinfectives	Hypnotics	Vaginal Products	Ophthalmic	
Antineoplastic Agents	Stimulants	Misc. GU Products	Otic	
Antineoplastics	Misc. Psychother.			

Clear User

Exit

Figure 16: Drug Class

(Image: drug-groups-er0141.bmp)

7. Click on the appropriate drug group button.  
(ex. Macrolide Antibiotics)

## Creating Custom Medication Menus for Administration in the ED (continued)

### Option 2: By Drug Classification(continued)

**Drug Picker** Active User: Thomas Daley

Patient Name:  Adm #  Pat #  SS #  Bed

Clindamycin HCl 150 MG (Clindamycin HCl Cap 150 MG)

Cleocin 300 MG (Clindamycin HCl Cap 300 MG)

Clindamycin HCl 300 MG (Clindamycin HCl Cap 300 MG)

Clindamycin HCl (Clindamycin HCl Powder)

Cleocin 75 M G/5ML (Clindamycin Palmitate For Susp 75 MG/5ML)

Clindamycin Palmitate 75 MG/5ML (Clindamycin Palmitate For Susp 75 Ml)

Cleocin 150 MG/ML (Clindamycin Phosphate Inj 150 MG/ML)

Clindamycin Phosphate 150 MG/ML (Clindamycin Phosphate Inj 150 MG/

Clindamycin Phosphate (Clindamycin Phosphate Powder)

Cleocin (300 ) (Clindamycin Phosphate in D5W IV Soln 300 MG/50ML)

Clindamycin Phosphate in D5W 300 MG/50ML (Clindamycin Phosphate in

Cleocin (600 ) (Clindamycin Phosphate in D5W IV Soln 600 MG/50ML)

Clindamycin Phosphate in D5W 600 MG/50ML (Clindamycin Phosphate in

Cleocin (900 ) (Clindamycin Phosphate in D5W IV Soln 900 MG/50ML)

Clindamycin Phosphate in D5W 900 MG/50ML (Clindamycin Phosphate in

Azithromycin 250 MG (Azithromycin Cap 250 MG)

**Zithromax 250 MG (Azithromycin Cap 250 MG)**

Biaxin 250 M G (Clarithromycin Tab 250 MG)

Clarithromyc in 250 MG (Clarithromycin Tab 250 MG)

Biaxin 500 M G (Clarithromycin Tab 500 MG)

Clarithromyc in 500 MG (Clarithromycin Tab 500 MG)

Biaxin 125 M G/5ML (Clarithromycin For Susp 125 MG/5ML)

Clarithromyc in 125 MG/5ML (Clarithromycin For Susp 125 MG/5ML)

Biaxin 250 M G/5ML (Clarithromycin For Susp 250 MG/5ML)

Clarithromyc in 250 MG/5ML (Clarithromycin For Susp 250 MG/5ML)

Drug Group:

Drug Subgroup:

Drug Class:

**Figure 17: Drug Picker**

(Image: drug-picker\_zit-er0143.bmp)

8. Select (*Double click*) the drug.

## Creating Custom Medication Menus for Administration in the ED (continued)

### Option 2: By Drug Classification (continued)

No.	Type	Route	Frequency	Duration	Name Displayed on Your Custom List
1/4	tab	po	1 time only	1 day	Zithromax 500 mg load dose
1/2	cap	pr	q day	2 days	Drug Name
3/4	mch	as	bid	3 days	Zithromax 250 MG
1	lozenge	od	tid	4 days	Generic
2	ampule	ou	qid	5 days	Azithromycin Cap 250 MG
3	packet	as	5 times daily	7 days	Drug Group
4	sup	ad	q 3*	10 days	Drug Subgroup
5	piece	au	q 3-4*	12 days	Drug Class
6	inplant	inhalations	q 3-6*	14 days	
7	patch	intranasal	q 3-12*	21 days	
8	bar	apply to affect. area	q 4*	28 days	
9	bottle	topically	q 4-6*	30 days	
10	gits	sublingual	q 4-8*		
<input type="checkbox"/> + 1/2	tsp	vaginal	q 5*		Quantity <input type="text" value="10"/> Days <input type="text" value="5"/> Remaining <input type="text" value="165"/>
	tbs	as irrigation	q 6*		SIG (Use Option Return to manually break lines)
	cc	transdermal	q 6-8*		2 cap po now
	mcg	IM	q 6-12*		
	mg	IV	q 8*		
	gm	subcutaneous	q 8-12*		
	ml		q 12*		
	mcu		<input type="checkbox"/> As Needed		
	mcu				
	ml				
	units				

To add 1/2 to the No. (i.e. 3-1/2) click the box on "+ 1/2" then click the appropriate No. button.

SIG may be built using buttons or typed.

The same drug can be repeated with different specifications. Use a different name for your custom list.

Save Change  
Save New  
Clear User  
Exit

Figure 18: Zithromax 500 Drug Screen

(Image: zithromax\_500-cr0149.bmp)

9. Edit the display name, drug name, generic, quantity, duration and SIG fields as required.

Save New

10. Click the Save New. (returns to custom drug list)

Exit

11. Click the Exit button.

## Modifying Custom Medication Menus for Administration in the ED

From the "Active Patient List" screen




1. Click the Utilities button.



2. Click on the Edit Doctor Drug Order Pref. List button.



3. Enter your password in the password prompt window  
4. Click on the Continue button.



**Figure 19: Custom Drugs (order menu)**

(Image: custom-drugs\_order-er0134.bmp)

5. Select (*click on*) the drug to be modified from the drug list at left of screen.



6. Click on the Change button.

## Modifying Custom Medication Menus for Administration in the ED (continued)

No.	Type	Route	Frequency	Duration	Name Displayed on Your Custom List
1/4	tab	po	1 time only	1 day	Zithromax 500 mg load dose
1/2	cap	po	q day	2 days	Drug Name
3/4	inch	os	bid	3 days	Zithromax 250 MG
1	lozenge	od	tid	4 days	Generic
2	ampule	ou	qid	5 days	Azithromycin Cap 250 MG
3	packet	as	5 times daily	7 days	Drug Group
4	sup	ad	q 3*	10 days	Drug Subgroup
5	piece	au	q 3-4*	12 days	Drug Class
6	implant	inhalations	q 3-6*	14 days	
7	patch	intranasal	q 3-12*	21 days	
8	bar	apply to affect. area	q 4*	28 days	
9	bottle	topically	q 4-6*	30 days	
10	qlts	sublingual	q 4-8*		
<input type="checkbox"/> + 1/2	tsp	vaginal	q 5*		Quantity <input type="text" value="10"/>
	tbs	as irrigation	q 6*		Duration <input type="text" value="5"/> Days Remaining <b>165</b>
	cc	transdermal	q 6-8*		SIG (Use Option Return to manually break lines)
	mcg	IM	q 6-12*		2 cap po now
	mg	IV	q 8*		
	gm	subcutaneous	q 8-12*		
	ml		q 12*		
	miu		<input type="checkbox"/> As Needed		
	au				
	iu				
	units				

To add 1/2 to the No. (i.e. 3-1/2) click the box on "+ 1/2" then click the appropriate No. button.

SIG may be built using buttons or typed.

The same drug can be repeated with different specifications. Use a different name for your custom list.

Save Change  
Save New  
Clear Uses  
Exit

**Figure 20: Zithromax 500 Drug Screen**

(Image: zithromax\_500-er0149.bmp)

9. Edit the display name, drug name, generic, quantity, duration and SIG fields as required.

Save Change

Exit

10. Click the Save Change button. (returns to custom drug list)

11. Click the Exit button.



## Deleting Custom Medications Menu Entries

From the "Active Patient List" screen



1. Click the Utilities button.



2. Click on the Edit Doctor Drug Order Pref. List button.



3. Enter your password in the password prompt window  
4. Click on the Continue button.





Figure 21: Custom Drugs (order menu)

(Image: custom-drugs\_order-er0134.bmp)

5. Select (*click on*) the drug to be deleted from the drug list at left of screen.



6. Click on the Delete button.



7. Click on the Exit button.

### **Tips and Hints: Medications Menu Entries**

- A different version of a drug can also be entered and Save New will save the additional version without changing the previous version.
- Note that 165 drugs can be listed in the custom drug list. The Remaining Available Positions notes how many positions remain.

## Managing Custom Prephrased Text

### Creating Custom Prephrased Text Items

From the "Active Patient List" screen

Utilities

Edit Prephrased Text

New Blank

1. Click the Utilities button.
2. Click on the Edit Prephrased Text button.
3. Enter your password in the password prompt window
4. Click on the Continue button. This will open the Physician Specific Prephrased Text.
5. Click the New Blank button at right of screen.
6. Select the appropriate medical record category from the pull-down menu.

Figure 22: Prephrased Text Editor (Physician Specific, new blank)

(Image: prephrased-text-editor\_1-er0366.bmp)

7. Enter desired text in the Title. Age and Text (center) fields as appropriate.
8. Click the Change button at right of screen.
9. Click on the Exit button at right of screen.

Change

Exit

## Modifying Custom Prephrased Text Items

From the "Active User" screen

**Utilities**

1. Click the Utilities button.

**Edit Prephrased Text**

2. Click on the Edit Prephrased Text button.

3. Enter your password in the password prompt window

4. Click on the Continue button. This will open the Physician Specific Prephrased Text.

**Prephrased Text Editor**

Last Edited By Welby      Date Last Edited 01/15/96

Title **Abdomen**      ☐ Male Only      ☒ Physician Specific  
 Min. Age:      ☐ Female Only      ☒ General Usage  
 Max Age:      **General Usage**

General Appearance  
 Mental Status  
 Skin  
 Head  
 Eyes  
 Ears  
 Mouth  
 Dental  
 Throat  
 Neck  
 Chest  
 Heart  
 Lungs  
**Abdomen**  
 Back  
 Genitourinary  
 Genitourinary  
 Rectal  
 Rectal  
 Extremities  
 Psych Admits  
 Neurological/Infant  
 Neurological/Child

**Physical Exam**

ABDOMEN: no tenderness, normal bowel sounds, no rebound, no guarding, no flank tenderness, no distention, not tympanitic. No peristaltic waves present. No ascites noted. No mass is palpable. No hernia is present. Liver and spleen are normal in size. Negative Murphy sign.

**Physician Specific**  
**General Usage**  
 New Blank  
 Change  
 Delete  
 Print All  
 Exit

**Figure 23: Prephrased Text Editor (General Usage, Abdomen selected)**

(Image: prephrased-text-editor\_2-er0366.bmp)

**Physician Specific**

or

**General Usage**

5. Click the Physician Specific or General Usage button at right of screen.

6. Select (**Double click**) the prephrased text title you wish to modify

7. Enter desired modifications to the text.

**Change**

8. Click the Change button at right of screen.

**Exit**

9. Click on the Exit button at right of screen.

### **Deleting Custom Prephrased Text Items**

*From the "Prephrased Text Editor" screen*



1. Select (*Double click*) the prephrased text title you wish to delete
2. Click on the Delete button.
3. Click on the Exit button at right of screen.

### **Tips and Hints: Prephrased Text**

- The General Usage list is editable and displayable by all physicians.
- The Physician Specific list is only editable and displayable by the physician account it was created under.


## TeleMed Data Maintenance

### **CAUTION**


*Do not attempt to use the Data Maintenance function if you are not experienced with it. Typically, this function will be used with the assistance of RLIS/TeleMed Support personnel.*

### **Accessing the Data Maintenance Screen**

*From the "Active Patient List"*



1. Click the Utilities button.
2. Click on the Maintenance button.
3. Enter your password in the password prompt window
4. Click on the Continue button.



**Figure 24: Maintenance**

(Image: maintenance-cr0184.bmp)

## Viewing Errors

*From the "Maintenance" screen*



1. Click the View Errors button.

The error log is a list of system errors to assist TeleMed programmers in locating problems. During phone conversations with TeleMed support personnel, the System Manager may be asked for information contained in the error log.

2. To exit, press your **Esc** key.

Date	Time	User	Error	Error

**Figure 25: Error Log Table**

(Image: error-log.bmp)

## Reindex (all) Data

*From the "Maintenance" screen*



**CAUTION:** Do NOT stop this process once started.

1. Click the Reindex All button.

Reindex All is a data recovery utility in the unlikely event that a data index has a problem. Although reindexing normally cannot damage any data, it should be used with care since, with a large accumulation of medical information, it may take a significant amount of time to complete. **Status** (refer to Figure 24: Maintenance, page 42) will indicate how many tables remain to be indexed and will indicate when the reindex process is complete.

## Testing TeleMed Changes

System managers can use the simulated patient function to enter pretend patients for teaching or experimentation. Entries made while in this mode will not be counted in any department statistics.

### Testing TeleMed with a Simulated Patient

*From the "Medical Information" screen*

**Physician**

1. Click the Physician button at the right of screen.

The screenshot shows the 'Patient ED Physician' interface. At the top, it displays 'Patient Name: Byrd, Newby' and 'Active User: J.E. Ross, MD'. Below this, there are fields for 'Adm #', 'Pat #', 'SS #', and 'Bed', all containing the value '01'. A 'Last Visit' field shows '1 / 1'. The main area contains an 'ED Physician' field with 'J.E. Ross, MD' and an 'A' button. To the right is a checkbox for 'Dictation Outstanding' and a 'Clear Physician' button. Below the physician field is a text instruction: 'The automatic button ("A") will assign the Active User as the ED Physician.' Underneath is an 'Assign ED Physician' label and a dropdown menu. At the bottom left is a checkbox for 'Simulated Patient'. At the bottom right are three buttons: 'Exit No-Save', 'Clear User', and 'Exit'.

**Figure 26: Patient ED Physician**

(Image: patient-ed-physician-er0233.bmp)

**Exit**

2. Check the **Simulated Patient** checkbox at lower left of screen
3. Click the Exit Button.



## Chapter 3: TeleMed Support Services

### TeleMed™ Service

Under the TeleMed Service program, RLIS provides support services to your local help desk and users. It's the practical, economical solution for hospitals with distributed information service groups and TeleMed software users. Three levels of support are offered under this program:

- **Platinum Key Support:** Includes toll-free help desk (standard and advanced), master copies of all maintenance releases, Quarterly Updates<sup>†</sup> and periodic upgrades<sup>†</sup> to the current version of your licensed software and On-Site<sup>‡</sup> Technical Support.
- **Gold Key Support:** Includes toll-free help desk (standard and advanced), Quarterly Updates<sup>†</sup> and periodic upgrades<sup>†</sup> to the current version of your licensed software and "fee per visit" on-site technical services.
- **Silver Key Support:** Provides all Quarterly Updates<sup>†</sup> and periodic upgrades<sup>†</sup> to the current version of your licensed software and "fee per call and per visit" technical services.

Support Provided	Platinum Key	Gold Key	Silver Key
<b>Standard Help Desk (toll free)</b> covering common operating errors, how-to questions and less complex problems.	24x7	8am - 5pm (CT)	Fee per call
<b>Advanced Help Desk (toll free)</b> For more complex problems, including those requiring program modification, RLIS senior staff is available between 8 AM and 5 PM Central time	8am - 5pm (CT)	8am - 5pm (CT)	Fee per call
<b>Remote Service (requires Remote Link Option)</b>	Optional	Optional	Optional
<b>On-Site Service</b> If necessary, RLIS will dispatch a technician to the customer's site	Yes	Fee per visit	Fee per visit
<b>TeleMed Maintenance Releases</b>	Yes	Yes	Yes
<b>TeleMed Quarterly Updates</b>	Yes	Yes	Yes
<b>TeleMed Periodic Upgrades</b>	Yes	Yes	Yes
<b>TeleMed Functional Upgrades</b>	No	No	No

- ♦ **Remote Link Discount Option:** When selecting the Remote Link Option, all sites that provide TeleMed dedicated ISDN connection(s) to their TeleMed servers will receive a discount on any of the three support options.

<sup>†</sup> see "Definitions" section for descriptions

<sup>‡</sup> as determined by TeleMed Service staff

### ***Contacting TeleMed Help Desk***

You can reach TeleMed support via phone or fax.

**Phone:** 800-496-7547

**Fax:** 210-495-8899

## Appendixes

### Appendix A: TeleMed License Agreement

RLIS and Licensee agree that the following terms and conditions will apply to all computer program products, user manuals, other documentation and services ("Program(s)") provided by RLIS to Licensee including Program(s) specified in any "Schedule" executed by both parties. The term Program(s) also includes any subsequent updates, modifications and enhancements to the Program(s) that are furnished by RLIS to Licensee.

1. License. RLIS hereby grants to Licensee a nonexclusive license (the "License") to use the Program(s) internally in the United States in accordance with the terms and conditions set forth in this Agreement and the relevant "Schedule." The License authorizes the use of the Program documentation furnished by RLIS and the use of any subsequent updates to the Program(s) that are furnished by RLIS to Licensee. The License authorizes the Licensee to install or use the copy of the Program(s) furnished to Licensee, and the back-up or archival copy made by Licensee, only on a single central processing unit ("CPU") (or authorized server(s) and workstation(s) under the network option license) identified by machine identification numbers provided by Licensee to RLIS. Installation or use of the Program(s) on more than one CPU at the same time requires a separate License for each such installation or use. Licensee may physically transfer the Program(s) from one CPU to another, provided that the Program(s) reside on only one CPU at a time. Licensee may not install the Program(s) on a network or electronically transfer the Program(s) from one CPU to another via a network, unless Licensee has obtained a network option for the Program from RLIS.

If RLIS delivers a network option for the Program(s) to Licensee, each such option shall consist of a network option license which shall authorize the Licensee to install and operate the Program(s) on the licensed server or workstation part of the Program(s), as the case may be, on a single authorized server or multiple workstations or clients identified by a machine identification number and located at a site provided by Licensee to RLIS so that the Program(s) may be accessed by users only in a network designated by Licensee only by means of the authorized server(s).

Unless otherwise specified in the relevant "Schedule", the maximum number of concurrent users for Program(s) licensed in a network option shall be limited to the maximum number of licenses granted to Licensee and in effect for such Program(s) or, in the case of updates to the Program(s) delivered under Support Services, the maximum number of licenses for which Licensee has purchased Support Services. Use of Programs in a network option by additional concurrent users shall require that Licensee purchase a commensurate number of additional licenses or, in the case of updates, Support Services for a commensurate number of licenses.

Licensee shall promptly notify RLIS of any changes in the machine identification number or location of the authorized server(s) or workstation(s) or CPU(s) from the one originally specified. If Licensee wishes to substitute a different authorized server or workstation or CPU under this Agreement, Licensee shall notify RLIS and RLIS shall either provide the Program(s) on the different authorized server or workstation or CPU or enable the existing Program(s) to operate on the different authorized server or workstation or CPU as soon as reasonably practicable and Licensee shall pay RLIS the amount, if any, specified in the relevant "Schedule" attached hereto.

The Licenses cannot be assigned, sublicensed or otherwise transferred to another person or legal entity without the prior written consent of RLIS, except that Licensee may permit its majority owned subsidiary to exercise any of Licensee's rights under this Agreement, provided that (i) such subsidiary agrees (in a writing provided to RLIS if so requested) and Licensee hereby represents that subsidiary will comply with all the terms and conditions of this Agreement, and it does so comply, (ii) Licensee shall be liable for any violations of the terms and conditions contained herein, and (iii) no such subsidiary is in either direct or indirect competition with RLIS as a provider of computer-aided software engineering products or products to develop, implement and maintain business critical information systems.

The Program(s) may contain an automatic shut-off or time out feature which will disable the Program(s) after a predetermined period of time if Licensee fails to provide RLIS with CPU, server or workstation or other machine identification numbers.

2. Term. This Agreement shall commence as of the date upon which this Agreement is signed by both parties and this Agreement and the License(s) shall continue unless and until terminated pursuant to the provisions set forth in this Agreement.

3. Charges. All fees are payable in US dollars. Upon receipt of the Program(s), the license and/or service fee payment terms set forth in the "Schedule" shall be due and payable. Payment shall be made at the office of RLIS, Inc. at 4319 Medical Drive #131-341, San Antonio, Texas 78229, or as otherwise directed by RLIS, Inc. in writing. Should Licensee request that delivery of the Program(s) occur more than thirty (30) days after execution of the "Schedule", such fees shall be subject to any changes in RLIS' commercial list prices. Licensee shall pay all sales, use, excise, and similar taxes applicable to the Program(s) and this Agreement and shall reimburse RLIS for any such taxes paid by RLIS, excluding taxes based on the net income of RLIS.

4. Delivery. RLIS shall deliver the computer programs to Licensee in the form of CD-ROM(s), diskette(s), or tape(s) containing machine readable code F.O.B. San Antonio, Texas, USA. In addition, RLIS shall provide Licensee with one set of RLIS' current user manual(s) for each computer program licensed hereunder describing the operation and use of that computer program.

5. Installation. RLIS shall install the computer programs as set forth in the Performance Criteria Checklist attached to this Agreement and executed by the parties.

6. Protection of Program. Licensee acknowledges that RLIS retains all title, copyrights, patents and other proprietary rights in and to the Program(s), including any and all copies, updates, modifications, enhancements, translations and other derivative works that duplicate or are based thereon. If RLIS provides consulting services, including on-site services or training services (collectively, the "Consulting Services") in connection with Licensee's use of the Program(s), such Consulting Services shall be provided at RLIS' then current rates. Licensee agrees any updates, modifications, enhancements, and translations to the Program(s) shall be the sole and exclusive property of RLIS, subject to use by the Licensee under the terms of this Agreement. Notwithstanding the foregoing, RLIS makes no claim of any proprietary rights in the output generated by Licensee's authorized use of the Program(s) or other code created by use of the Program(s).

Licensee shall not make copies of the Program(s), except Licensee may make a back-up or archival copy of each computer program and such additional copies as may be required by law. Additionally, if a Program is delivered in CD-ROM form, Licensee may, for Licensee's internal use only, reproduce the Program object code on diskette or tape or reproduce the user manual in printed or other form, provided that the maximum number of concurrent users of such Program, and the maximum number of copies of such user manual made by Licensee, are limited in each case to the maximum number of licenses granted to Licensee and in effect for such Program or, in the case of updates to the Program delivered under Support Services, the maximum number of licenses for which Licensee has purchased Support Services. Licensee shall not remove any copyright or other proprietary rights notice included in or on any Program or user manual, and shall reproduce all such notices in or on all copies made by Licensee.

Licensee shall limit access to the Program(s) to its employees and agents whose responsibilities require such access, and Licensee shall adopt reasonable measures to assure that its employees and agents will make no disclosure of the Program(s) or the information contained therein to other persons or legal entities. Licensee agrees to treat any computer program code and documentation furnished to Licensee as a confidential trade secret and valuable asset of RLIS and agrees that such code shall not be used for any purpose other than to assist in the normal use of the Program(s) as defined in the user manuals provided with the Program(s). In particular, but without limitation, Licensee agrees it will not decompile, disassemble or attempt in any way to reverse engineer or access the source code, excluding such access to source code as may be required by law, or to develop a competing product based on the Program(s). Licensee shall not rent, loan, or permit any timesharing or service bureau use with respect to the Program(s).

7. Warranties, Disclaimer and Limitation. RLIS warrants that it has the right to grant the License(s) and that the Program(s) furnished to Licensee does not infringe any United States patent, copyright or trade secret rights of any third party. RLIS shall, at RLIS' expense, defend or settle any action brought against Licensee based on a claim that the Program(s) infringe a United States patent, copyright or trade secret right and will pay all costs and damages finally awarded against Licensee in any such action which are attributable to such claim; provided that (i) Licensee promptly notifies RLIS in writing of the claim; (ii) RLIS shall have sole control of the settlement and defense of any action to which this indemnity relates; and (iii) Licensee cooperates with RLIS in every reasonable way to facilitate such defense or settlement. If a Program becomes or, in opinion of RLIS, is likely to become the subject of a claim of infringement, RLIS may, at its option, procure for Licensee the right to continue using the Program, replace or modify it to make it non-infringing or, if it determines that neither of the foregoing alternatives is reasonably available to RLIS, then RLIS may terminate the License(s) for the Program upon forty-five (45) days notice to Licensee and refund the pro-rated (based on a 3 year depreciation "Schedule") license fee paid for the Program, and, except as set forth in this paragraph, RLIS shall leave no further liability to Licensee. RLIS shall not be liable for any infringement or claim thereof based upon the use of the Program(s) or any element of the Program(s) in combination with computer programs or hardware not supplied by RLIS, or upon any modification to the Program(s) made by any party other than RLIS.

THE FOREGOING PARAGRAPH OF THIS SECTION STATES THE ENTIRE LIABILITY OF RLIS FOR ANY LOSS AND DAMAGES WHATSOEVER AS A RESULT OF THE INFRINGEMENT OF ANY COPYRIGHT, PATENT, TRADE SECRET OR OTHER INTELLECTUAL PROPERTY RIGHTS.

RLIS warrants that each copy of the Program provided by RLIS will be free of physical defects in the media when delivered to Licensee. RLIS shall replace any copy that does not comply with this warranty, provided that Licensee notifies RLIS of the defect within thirty (30) days of the delivery date. Replacement shall be Licensee's sole remedy for any defective media.

THE EXPRESS WARRANTIES SET FORTH IN THIS SECTION ARE IN LIEU OF ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED, INCLUDING THE WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE.

Licensee is solely responsible for: the selection of the Program(s) and the appropriate computer equipment for Licensee's needs. The proper use of the Program(s) by qualified personnel, taking measures to prevent loss of data and all other matters under its control. Except with regard to the warranty provided in the first paragraph of this section 6, RLIS' liability and contract, tort or otherwise with respect to any Program(s) and/or Support Services, shall not exceed the amount of the license fee(s) and support fee(s) paid to RLIS by Licensee with respect to such Program(s) during the preceding twelve (12) months.

IN NO EVENT SHALL RLIS BE LIABLE FOR SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES, INCLUDING ANY DAMAGES RESULTING FROM LOSS OF USE, LOSS OF DATA, LOSS OF PROFITS OR LOSS OF BUSINESS, EVEN IF RLIS HAS BEEN ADVISED OF THE POSSIBILITY OF SUCH DAMAGES.

8. Support Services. RLIS will provide the Support Services (as hereinafter defined) for the Program(s) free of charge during the Warranty Period, if any, set forth in the "Schedule" applicable to such Program(s) and at the then-current Support Service fee(s) thereafter. Prior to the expiration of the Warranty Period, if any, and at each Support Service renewal thereafter, Licensee may select one of the Support Services options in RLIS' then current published price list ("TeleMed Services") and RLIS shall provide the Support

Services at the then current fee(s) in acceptance with the terms of the selected option as set forth in such price list. Updates to the Program(s) shall not extend the term of any Warranty Period.

**WARNING:** Failure to implement Quarterly Updates within the 45 day time period will result in automatic shutdown of the prescription and medications components of the TeleMed system.

Licensee assumes all risks arising from Licensee's failure to implement updates and any other corrections released by RLIS and from implementation of Program modifications not supplied by RLIS. Support Services will be provided by mail, telephone or on-site as specified under the Support Services option selected by Licensee.

Support Services for the first year of Support Services purchased by Licensee following the Warranty Period will be due forty-five (45) days prior to closure of Warranty Period. Support Services are renewable on an annual basis thereafter at RLIS' then-current applicable Support Services fees, which will be due on or before each anniversary of the expiration of the Warranty Period. If Licensee fails to renew Support Services, RLIS shall thereafter have no obligation to correct or support the Program(s). Notwithstanding the above, RLIS may at any time, upon one year's notice to Licensee, terminate Support Services for any Program(s).

9. **Quarterly Updates.** Quarterly Program updates are distributed by RLIS to implement changes in prescription, standard forms or medical information. These updates will be distributed the first business day of the following months: January, April, July, October. Quarterly Updates must be installed within 45 days of distribution.

**WARNING:** Failure to implement Quarterly Updates within the 45 day time period will result in automatic shutdown of the prescription and medications components of the TeleMed system.

**NOTE:** In the interest of practicing good medicine, these Quarterly Updates address: pharmaceutical database updates, 3<sup>rd</sup> party licensing terms, medical knowledge base updates and other enhancements as required for prudent medical practice.

10. **Termination of License.** RLIS may terminate this Agreement or the License for any and all Program(s), effective immediately upon written notice to Licensee, if Licensee fails to perform any of its obligations under this Agreement and fails to remedy such breach within thirty (30) days after RLIS gives Licensee written notice of the breach.

In the event that this Agreement for a License is terminated for any reason, Licensee shall promptly return to RLIS the original and all copies of the Program(s), including all documentation covered by the terminated License(s). The provisions set forth in section 5, 6, 8 and 12 of this Agreement shall survive any termination of the License.

11. **U.S. Government Restricted Rights.** If Licensee is a U.S. government agency or government subcontractor, Licensee hereby acknowledges that the computer programs comprising the Program(s) are commercial computer software developed at private expense and provided to Licensee subject to RESTRICTED RIGHTS, and all other Programs' technical data, including without limitation any user manuals and other documentation, are provided to Licensee subject to LIMITED RIGHTS. The terms "commercial computer software," "private expense," "RESTRICTED RIGHTS" and "LIMITED RIGHTS" as used in this section 9 are defined in 48 CFR: DFAR 252.227-7013, FAR 52.227-19 and DFAR 52.227-7013.

Licensee agrees that it will not remove or permit to be removed the following legend and notice printed on the Program(s) and will place this legend and notice on all copies made by Licensee:

RESTRICTED RIGHT LEGEND - Use, duplication or disclosure by the Government is subject to restrictions as set forth in subparagraph (c) (i) (ii) of the Rights in Technical Data and Computer Software clause at DFAR 252.227-7013 (48 CFR). Contractor/manufacturer is RLIS, Inc., 4319 Medical Drive #131-341, San Antonio, Texas 78229

NOTICE - Notwithstanding any other lease or license agreement that may pertain to, or accompany the delivery of this restricted computer software, the rights of the Government regarding use, reproduction and disclosure are as set forth in subparagraph (c) (1) and (2) of Commercial Computer Software - Restricted Rights clause at FAR 52.227-19.

12. **Exportation.** Licensee shall comply with all applicable laws, including, without limitation, the export control laws of the United States of America and other countries and the relevant regulations issued by the United States Department of Commerce and Department of State, concerning the import, export and re-export of the Program(s), any other accompanying documentation and the direct products thereof. Without limiting the generality of the foregoing, Licensee shall not, and Licensee hereby assures RLIS that it will not, allow the export or re-export, directly or indirectly, of any Program, accompanying documentation or the direct products thereof unless prior written authorization is obtained from RLIS and, where required, the United States government.

13. **Third Party Products.** If the Program(s) include any third party computer programs or other products, then the following terms shall apply to such third party products and in the event of a conflict between this Agreement and the following terms, the following terms shall govern with respect to third party product: (i) any shrink-wrap agreement or other terms or conditions included in the third party product packaging or documentation; and (ii) any terms or conditions applicable to third party products that are attached to, or incorporated into, a "Schedule."

14. **Severability.** If any provision of this Agreement, or portion thereof, or the application thereof to any circumstance shall be held to be invalid or unenforceable, the remainder of this Agreement and the application thereof to other circumstances shall be nevertheless valid. In lieu of such invalid or unenforceable provision, there shall be added automatically a provision as similar in terms to such invalid or unenforceable provision as may be possible and be legal, valid and enforceable

15. General. This Agreement sets forth the entire agreement and understanding of the parties relating to the subject matter hereof and supersedes any and all oral and prior written agreements, understandings and quotations relating thereto. No alteration, modification or cancellation of any of the provisions of this Agreement shall be binding unless made in writing and signed by officers of the parties. Printed terms and conditions on Licensee's Purchase Order(s) shall not apply to Program(s) obtained hereunder. This Agreement will be governed by, and construed and enforced in accordance with, the substantive law of the State of Texas, USA. The UN Convention on the International Sale of Goods shall not apply. The English language versions of this Agreement shall govern.

This Agreement shall be binding upon and inure to the benefit of the parties and their respective successors, permitted assigns and legal representatives.

# Index

## A

Account  
   deleting, 16  
   modifying, 15  
 Account Management, 14, 15

## B

Backup, 9

## C

Creating  
   custom medications menus, 29  
   custom medications menus, by  
     drug classification, 32  
   custom medications menus, by  
     drug name, 30  
   custom prephrased text, 39  
   custom prescription menus, 19  
   custom prescription menus, by  
     drug classification, 23  
   custom prescription menus, by  
     drug name, 20  
   medical staff directories, 17  
   phone directories, 17  
 Custom Medications, 29  
 Custom Prescription Menus, 19

## D

Data Links  
   to other hospital systems, 8  
 Defragmentation, 10  
 Deleting  
   accounts, 16  
   custom medications menus, 37  
   custom prephrased text, 41  
   custom prescription menus, 27  
   medical staff directories, 18  
   phone directories, 18

## E

Edit

custom medications menu  
   entries, 35  
 custom prescription menu  
   entries, 25, 26, 36  
 phone directory entries, 18  
 prephrased text, 39  
 Errors  
   viewing, 43

## H

Help Desk, RLIS/TeleMed, 46

## I

Installation Specification, 7, 8, 9,  
 10, 11

## L

License Agreement, 47

## M

Maintenance  
   data maintenance screen, 42  
   reindex data, 43  
   testing changes, 44  
   viewing errors, 42, 43  
 Managing the TeleMed System, 9  
   backup, 9  
   system startup, 9  
   updates & upgrades, 9  
 Modifying  
   accounts, 15  
   custom medications menus, 35  
   custom prephrased text, 40  
   custom prescription menus, 25,  
     26, 36  
   medical staff directories, 18  
   phone directories, 18

## P

Password  
   modifying, 15  
 Phone Directories, 17

creating, 17  
 deleting, 18  
 modifying, 18  
 Prephrased Text, 39  
 Printing, 9  
   printer configuration, 9  
 Privileges, 16

## R

Reindex data, 43

## S

Security  
   overview, 7  
   TeleMed data, 8  
   TeleMed system application, 7  
 Simulated patient, 44  
 Storage and Archive Locations, 8  
 Support Services, 45  
 System Manager  
   rights, 15  
   Rights, 16  
 System Performance, 10  
   defragmentation, 10  
 System Startup, 9

## T

Tasks, typical, 13  
 TeleMed Data System Features, 8  
 Telephone Access, 7  
 Testing  
   changes, 44

## U

Updates & Upgrades, 9  
 User Rights, 14, 15







15600 San Pedro  
Suite 203  
San Antonio, Texas 78232  
Phone: 210-490-1800  
Fax: (210) 495-8899

**THIS PAGE BLANK (USPTO)**